

BEHANDLUNGSSTUHL FÜR PATIENTEN – HYD



Bestimmungsgemäße Verwendung

Der Patientenbehandlungsstuhl ist so konzipiert, dass er den Patienten in einer sitzenden oder zurückgelehnten Haltung stützt und positioniert, um einen einfachen Zugang und Patientenkomfort während eines täglichen klinischen Eingriffs wie Blutspende, Impfung, postoperative Behandlung und andere zu gewährleisten. Mit Zubehör kann der Stuhl für verschiedene Behandlungen verwendet werden.

- 1 Rückenteil**
- 2 Sitzteil**
- 3 Beinteil**
- 4 Doppelrollen Ø100mm mit Bremse**
- 5 Hydraulische Säule**
- 6 Linke Armlehne**
- 7 Rechte Armlehne**

Optional

- 1 Kopfstütze
- 2 Papierrollenhalter

Specifications:

Technical data

Product code:	11000410
Length:	830 mm
Width:	610 mm
Width with armrests:	830 mm
Height:	1320 – 1600 mm
Lying surface length:	1820 mm
Seat section width:	580 mm
Seat section height:	570 - 850 mm
Back section length:	880 mm
Seat section length:	460 mm
Leg section length:	420 mm
Back section inclination:	-10° to 73°
Fixed seat inclination:	7°
Leg section inclination:	-90° to 7°
Trendelenburg:	-10°
Double-wheel casters:	Ø100
Weight:	45 kg
Safe working load:	170 kg
Armrest safe working load:	10 kg







BY AGREEMENT

The minimum order quantity is 5 pieces.

